

1. Technical Specification

Bed Data

| | |
|---|------------------------------|
| Overall Length without scales | 2270mm |
| <i>Extended</i> | 2355mm |
| Overall Length with scales | 2425mm |
| <i>Extended</i> | 2510mm |
| Overall Width | 1105mm |
| <i>Extended</i> | 1405mm |
| Mattress Platform Length | 2015mm |
| <i>Extended</i> | 2225mm |
| Mattress Platform Width | 915mm |
| <i>Extended</i> | 1220mm |
| Mattress Platform Height | |
| <i>Without scales</i> | Min 415mm - Max 835mm |
| <i>With scales</i> | Min 415mm - Max 835mm |
| Backrest – Max tilt angle | 55° |
| Upper Leg Section – Max tilt angle | 25° |
| Lower Leg Section – Max tilt angle | 18° |
| Trendelenburg | 21° |
| Reverse Trendelenburg | 21° |
| Safety Rail Height Above Mattress Platform | 460mm |
| Product weight without scales (approx) | 190kg |
| Product weight with scales (approx) | 215kg |
| Maximum Safe Working Load (SWL) | 413kg / 65 stone / 910lbs |
| Medical Device Classification – Without Scales | Class 1 |
| Medical Device Classification – With Scales | Class 1 Measuring |

Electrical Data

| | |
|---------------------------|---------------------------------|
| Power Supply | 230V~ 50-60Hz |
| Power Consumption | Max. 1.6A AC |
| Duty Rating | Intermittent 10% (2min / 18min) |
| Equipment Classification | Class II Type B Applied Part |
| Liquid Ingress Protection | IPX6 |