

bariatric commode

standard bariatric commode



Key features / spec

337kg (53 st) swl

710mm (28") seat width

Includes removable waste container with lid

Swing away arms enables side transfers

Large space between seat & rear handle for wider patients

Plastic bed pan insert (optional)

Our standard bariatric commode can also be used as a shower chair and will fit over a standard sized toilet, making the commode extremely versatile.

Made with an aluminium frame, it has four strong braked castors to ensure stability whilst in use and the swing away arms enable easy side transfers.

Product code: -
BMC002

Emergency 24/7 Rentals

Supporting Bariatric Care



t: 0333 800 9000

f: 0333 800 9001

e: info@benmormedical.co.uk

w: www.benmormedical.co.uk