

bariatric drop end couch




Low patient transfer height

Width
105cm



38cm 49cm 93cm



Height range 50 to 105cm



Key Features

Safe working load and lifting capacity of
325kg (51st)

Memory foam upholstery for additional
patient comfort

Adjustable angle electric backrest
operation from 0° to +80°

Electric height elevation from 50cm to
105cm

Electric foot sections, adjustable in angle
from 0° to -90°

Purpose made precision bearings on all
pivot points

The bariatric drop end couch feature electric independently operated foot sections for effortless positioning for various lower limb procedures. The electric backrest supports the patient in the sitting position. There is a choice of models with either single or dual foot sections, all with a safe working load and lifting capacity of 325Kg (715lbs).



Optional Accessories

- Breathing hole and plug for prone treatments
- Head support cushion
- Side support rails (Pair)
- Side support rail cushions (Pair)
- Paper roll holder