

bariatric wheelchair

folding bariatric wheelchair



Key features / spec

210kg, 254kg and
320kg swl available

Folds for easy storage
and transportation

Removable, adjustable footplates

Detachable arms and
padded arm rests

Removable seat
and back cushions

Built in strengthening and bracing to
cope with heavier patients

Our range of folding wheelchairs consists of a variety of transit and self-propelled wheelchairs which are available with seat widths and depths to accommodate most patients.

Front and back seat supports use telescopic guides to ensure the seated area is guided to its supports. Built in strengthening and, in particular, bracing to reduce the possibility of twisting of the frame with a minimum of moving parts to ensure a user friendly wheelchair.

Product codes: -

BMWC001 - 210kg (33st) self-propelled

BMWC002 - 210kg (33st) transit

BMWC003 - 254kg (40st) self-propelled

BMWC004 - 320kg (50st) self-propelled

Emergency 24/7 Rentals

Supporting Bariatric Care



t: 0333 800 9000

f: 0333 800 9001

e: info@benmormedical.co.uk

w: www.benmormedical.co.uk