

# bariatric rollator

adjustable bariatric rollator



## Key features / spec

320kg (50st) swl

Strong steel construction with a durable and reliable frame

Folds for storage and easy transportation

Width adjustable from 635mm to 760mm (25" - 30")

Height adjustable from 510mm - 915mm (20" - 36")

Pressure activated braking system

Product code: -  
BMR001

Our bariatric rollator offers stability, safety and freedom of movement for the larger patient and can meet the needs of people who have difficulty with their balance.

The rollator is constructed from steel with strong welded joints. It is lightweight and easy to store. The castors are skid resistant and easy rolling with a braking system on the rear wheels which is activated by applying weight to the frame.

A range of accessories are also available including detachable seats, baskets and gutter supports for the patient's arms.

Emergency 24/7 Rentals

Supporting Bariatric Care



t: 0333 800 9000

f: 0333 800 9001

e: [info@benmormedical.co.uk](mailto:info@benmormedical.co.uk)

w: [www.benmormedical.co.uk](http://www.benmormedical.co.uk)