

cathedra static chair

bariatric static chair



Key features / spec

Available with a 254kg and 320kg
(40st and 50st) swl

Available with 610mm, 710mm and
815mm (24", 28", 32") seat widths

Breathable and waterproof
fabrics for easy cleaning

Rear castors for easy
transportation and position

Large choice of fabrics

Seat height adjusters, easy release
drop and padded arms available

We offer a selection of static hospital chairs which are available in a variety of finishes, weight limits and seat widths.

Reflexion foam cushions are fitted as standard for medium risk pressure relief. Chairs are supplied with housekeeping wheels as standard to enable easy transportation. Optional extras include drop arms for easy side transfers, seat height adjusters and arm pads for added patient comfort.

Product codes: -

BMSC001 - 610mm (24") seat width, 254kg (40st) swl

BMSC002 - 710mm (28") seat width, 320kg (50st) swl

BMSC003 - 815mm (32") seat width, 320kg (50st) swl

Emergency 24/7 Rentals

Supporting Bariatric Care



t: 0333 800 9000

f: 0333 800 9001

e: info@benmormedical.co.uk

w: www.benmormedical.co.uk