

# deporto

bariatric commode and transfer chair



- Wheeled frame
- Load capacity 400kg (60st)
- 28" or 31" (711mm or 787mm) seat width
- 17" (431mm) seat height
- Powder coated steel
- Padded seat and back
- Flip up armrests
- Removable foot complete with footplate
- Strong footplate that is used to support the patient
- Removable seat and backrest for easy cleaning-infection control compliant

features	benefits
400kg (60st) swl	Designed to support heavy loads
28" or 31" (711mm or 787mm) seat width	Wide enough for bariatric patients
Adjustable 'pram handles'	Ensure carer best positioning when pushing patient
Strong footplate	To provide support for the patient
Padded seat and back	For extra comfort for the patient
17" (413mm) seat height	Can fit over most standard toilet seats