

dynamic cushion

pressure relieving air cushion



Key features / spec

Inboard zip, minimises risk of contamination, avoids ingress of fluids

Two way stretch material, conforms to the contours of the patient's body

Various widths available, compatible with our seating range

Adjustable pressure, provides optimum pressure relief for patients

72 hours air retention when pump is disconnected

Waterproof, wipe clean cover for easy cleaning and infection control

Our dynamic air cushion has been designed to give high risk pressure relieving care to bariatric patients who suffer from pressure sores and general discomfort.

The cushions have an alternation cycle of ABAB and operate using a pump with adjustable pressure. The cover is made from a two way stretch material which conforms to the contours of the body and is vapour permeable, water resistant, fire retardant and soft whilst still being durable.

Product codes: -

BMAC001P - 286kg (45st) swl

BMAC002P - 349kg (55st) swl

BMAC003P - 349kg (55st) swl

Emergency 24/7 Rentals

Supporting Bariatric Care



t: 0333 800 9000

f: 0333 800 9001

e: info@benmormedical.co.uk

w: www.benmormedical.co.uk