

# portable bed scales

600kg weight capacity

The bed scales has a capacity to weigh patients up to 600kg and includes a tare function to deduct the weight of the bed.

The lightweight beams are simple to use and provide a straightforward weighing solution for bed dependant patients. The two lightweight beams have handles and wheels for ease of transportation.



- Easy and accurate weighing of bedridden patients
- Portable bed weighing beams
- Tare to deduct weight of bed



Capacity	600kg
Graduations	200g
Dimensions	800mm x 120mm x 15mm
Weight of scale	11.5kg
Power Supply	Mains / Rechargeable



Product code: BMW004