


portable floor scale

extra wide low profile platform with 500kg weight capacity

The portable floor scale allows for discrete weighing of larger patients that traditional scales found in GP surgeries cannot.

The extra wide base enables obese and elderly patients, unable to stand on a small platform, to be weighed safely and accurately. The lightweight aluminium construction and the carry case make the scale ideal for home visits.



- Extra wide, low profile suitable for all patients
- Discreet weighing
- High capacity
- Ideal for GP surgeries & weight management clinics
- Light platform

Easily stored in carry case

Capacity	500kg
Graduations	100g<300kg>200g
Dimensions	600mm x 50mm x 390mm
Weight of scale	8kg
Power Supply	Mains / Rechargeable



Product code: BMW002