

rota stand compact

bariatric rota stand



Key features / spec

200kg (32st) swl

Foldable and compact, reduces storage space needed

Folded length 1140mm, folded width 525mm

Only weighs 10kg (2st)

Height adjustable frame and shin pad to suit the majority of patients

Stabilizing bar allows a one person transfer

The rota stand compact combines all the benefits of a standard rota stand but with the added value of being foldable yet still lightweight.

The nature of the stand means that it is ideal for community assessments by therapists who often have to lift equipment in and out of their vehicles. Equally, a more independent lifestyle is encouraged for clients who wish to travel to friends and relatives.

The rota stand compact fits easily into the majority of car boots.

Product codes: -
BMHE051C

Emergency 24/7 Rentals

Supporting Bariatric Care



t: 0333 800 9000

f: 0333 800 9001

e: info@benmormedical.co.uk

w: www.benmormedical.co.uk