

stand on scale

suitable for bariatric patients up to 500kg

The stand on scale provides a large weighing platform and low profile access suitable for bariatric and geriatric patients.

This scale is ideal for hospital wards, outpatient departments and weight management clinics. It has a large capacity of 500kg so it is useful for weighing bariatric patients. The handrails make it easy to weigh patients that are unable to stand unassisted.



- Curved, strong handrails for supporting patients
- Suitable for bariatric patients
- Low-level platform for easy access
- Wheels for easy movement



Capacity	500kg
Graduations	200g
Dimensions	700mm x 600mm x 50mm
Weight of scale	33kg
Power Supply	Mains / Rechargeable



Product code: BMW006