

verto turning mattress

bariatric turning mattress



Key features / spec

455kg (71st) swl

2032mm x 1220mm x 254mm
(L x W x H)

Mattress weight 11kg

Automatic lateral rotation, reduces
the risk of nursing staff injury

Adjustment pressure settings, enhances
patient comfort and pressure relief

Easily accessible lock key,
increases patient safety

Product code: -
BMPB071TM

Many patients suffering from the ability to move or turn themselves develop pressure sores where insufficient blood is supplied to various areas of the body.

The gentle, gradual, 30-40° turning angle of the patient from side to side by the verto turning mattress exerts less than 26mm hg of pressure on the skin tissue.

Air escapes through tiny pinholes in the top of the cells to circulate under the cover. The moisture from the patient transfers through the vapour permeable fabric cover to interact with the air blowing up from the mattress.

Emergency 24/7 Rentals

Supporting Bariatric Care



t: 0333 800 9000

f: 0333 800 9001

e: info@benmormedical.co.uk

w: www.benmormedical.co.uk